



# Hot Plate Safety

**Hot plates used in labs present many potential dangers, such as fires and electrical shock. Fires can cause injuries, significant disruption of lab operations and loss of scientific data.**

## Basic Precautions:

- Periodically test the function of the “off” switch on each hot plate to verify that it works and the device cools quickly when the switch is in the “off” position.
- Always inspect equipment prior to use. Do not use if the plug or cord is worn, frayed, or damaged, if the grounding pin has been removed, or if a spark is observed. Check for corrosion of the thermostat, which can also cause a spark.
- Do not use hot plates in the presence of flammable or combustible materials. Fire or explosion may result.
- Provide secondary containment for any flammable liquids being heated to prevent liquid from contacting the hot plate in the event of a leak or overflow.
- Read the manufacturer’s instructions before using, and register the device with the manufacturer so you will be notified of any warnings or recalls.
- **TURN OFF** the hot plate when not in use. The surface of a hot plate stays hot for some time – and looks exactly the same as a “cold” plate.

## Tips for Working Safely with Hot Plates:

- Use only heat-resistant, borosilicate glassware, and check for cracks before using. Do not place thick-walled glassware, plastic containers, soft-glass bottles or jars on a hot plate.
- The hot plate surface should be larger than the vessel being heated.
- Ensure that electrical cords and temperature sensor probe wires do not come in contact with the hot plate surface.
- Use a medium to medium-high setting to heat most liquids, including water. Do not use the high setting to heat low-boiling liquids. The hot plate surface temperature can reach up to 540°C.
- Do not place metal foil or metal containers on the hot plate - the top can be damaged and shock hazard may result.
- Use thermal gloves or tongs when removing hot items from a hot plate and allow items to cool prior to handling.

